



PtoleMedic System

**Instructions for Use - Uploading MRI to
Database (MRI Center)**

Software Release 1.3.x

Document ID: LB-72-01-002

Model: Revision A

January 24, 2020

Author:

Signature on File

Signature: _____

_____ **Date**

Reviewed and Approved by:

Signature on File

Signature: _____

_____ **Date**

This is the reference guide for the MRI Technologist as an aid in uploading the MRI scanning images to the Lento Medical Innovation cloud database. There is a separate full MRI reference guide for use in positioning and patient orientation. Also, a quick MRI setting reference guide is available for use by experienced technologists and both are available upon request if they have not already been supplied or is misplaced.



WARNING: Please note that MRI scanner in patients with metallic implants in or near the knee joint may adversely affect the quality and accuracy of the images obtained. It is recommended that MRI scans not be attempted.



WARNING: The use of custom cutting guides in pediatric patients has not been studied and the results of the use of this product in these patients is unknown. We do not recommend that such surgery not be attempted.



Manufacturer:

Lento Medical Innovation, Inc.
15110 Northwest Freeway, Suite 150
Houston, TX 77040 USA
+1 (510) 413-3230



The following organization is the Authorized Representative for the PtoleMedic System:



HealthLink Europe BV
De Tweeling 20-22
5215 MC's - Herrogenbosch
The Netherlands
Tel.: +31-13 5479300
Web: www.healthlinkeurope.com

Introduction

These instructions are intended to guide the MRI Center to upload MRI files through the Lento Medical Innovation web portal (www.lentomedical.net).

1. Log In

Lento Medical Admin or Sales personnel will assist the MRI representative with registering an account for the Lento Medical Innovation web portal. Once a username and password are created, the surgeon will use the information to log in.

- a. Enter Username
- b. Enter Password
- c. Click "Login" button

ProleMedic Surgical Planning Software Ver. 1.3.0

Please login to see this page.

Username:

Password:

[New User / Reset Password](#)



2. MRI Center Landing Page

User landing page will appear listing all surgeons, see below.

Surgeons

Doctor's Name --	Doctor's Office Name --	Doctor's Office Address --
Lento, Surgeon	Surgeon Lento Office	401 Kato Terrace, Fremont, CA 94539, USA

3. MRI Uploading Page

If you click name of surgeon, it will show all active cases for user to view and upload. The status of a case show "Ready to Upload MRI", see below

Case # --	Patient name	Date Of PO --	Date Of Surger...	Date Of Birth	Surgery type	Surgery side --	Status --
K58874R	Lento, Patient 29	11/04/2019	01/01/2020	09/13/1988	TKA	Right	<input type="button" value="Ready to Upload M"/>

User can click the case number and the MRI upload page will appear, see below.

K58874R Ready to Upload MRI

Dr:	Surgeon Lento
Office Name:	Surgeon Lento Office
Office Address:	401 Kato Terrace, Fremont, CA 94539, USA
eMail:	s1@lentomedical.com
Office number:	5104137014
Fax number:	5104137020

Patient Information	
Name:	Patient 29 Lento
Sex:	Female
Age:	31
Date Of PO:	11/04/19
Date Of Birth:	09/13/88
Date Of Surgery:	01/01/20
Surgery Type:	Total Knee Arthroplasty (TKA)

Drop Zone

Note: before uploading the MRI file, put all the MRI files “Coronal Kee”, “Coronal Hip”, and “Coronal Ankle” into one folder.

4. Uploading MRI

There are two ways you can upload MRI files.

A. Drag and drop

You can drop the MRI file folder into “Drop Zone”

B. Manual folder selection

You can click “Drop Zone” to select the folder and click “Upload” button.

5. Uploading Progress

You will see uploading progress bar as shown below.

Lento Medical Inc. [LOG OUT](#)


K58874R Ready to Upload MRI

Dr:	Surgeon Lento
Office Name:	Surgeon Lento Office
Office Address:	401 Kato Terrace, Fremont, CA 94539, USA
eMail:	s1@lentomedical.com
Office number:	5104137014
Fax number:	5104137020


51%

6. MRI upload confirmation

Once MRI files are uploaded correct, confirmation popup as shown below.




Once you click “OK”, you will be re-directiing to case list page as shown below.:



Note: once MRI is uploaded successfully, case no longer show in above list.

7. Log Out

From top menu, click “Log Out” menu and clock “Log Out” to log out from system, see below.



8. Follow up

Our technican will review the MRI files that had being upload and contact you if there is any issues.