PtoleMedic System

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
Slice Plane:	Coronal
Slice thickness (mm):	4
Spacing/Gap (mm):	1
Number of slices:	14
FOV(mm):	240 mm
Matrix:	256 x 160
NEX:	2
NPW/Anti-Aliasing/Fold Over Suppression:	ON
File Series Name:	"Coronal Hip"

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

The image shows the Lento Medical login page. At the top is the Lento Medical logo with the tagline 'MEDICAL INNOVATION'. Below the logo, it says 'PtoleMedic Surgical Planning Software Ver. 1.3.1'. The main text says 'Please login to see this page.' followed by input fields for 'Enter Username' and 'Enter Password'. There is a link for 'New User / Reset Password' and a 'LOG IN' button. At the bottom, there is a CE mark and the text '* Not available for sale in the U.S.A.'.

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The image shows the Lento Medical case creation form. At the top, it displays 'Case ID: K4fen6L' and 'Status: Waiting for MRI'. There is a 'COMPLETE' button. The form is divided into several sections: Patient, Surgeon, Distributor, Hospital, and MRI Center. Each section contains various input fields for personal and professional information. The Patient section includes fields for First Name, Middle Name, Last Name, Date of Birth, Sex, Age, Date of Surgery, Surgery Type, Surgery Side, Default Alignment, Date of Purchase Order, Implant Information, Ship To, Bill To, and July Distributor. The Surgeon, Distributor, Hospital, and MRI Center sections each have fields for Name, Company Name, Contact Person, Address, Email Address, Office Number, Mobile Number, and Fax. At the bottom, there is a 'COLLAPSE' button and a 'SURGEON APPROVED' button.

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HOME LOGIN

Surgeon Orders

User: Aug, July
Office Name: July MRI Center
Office Address: Third Street
Berlin
7454654, DE
Email Address: woo4@lento-medical.com
Phone Number: 49132132132
Fax Number: 4

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Surgeon List

Last Name	First Name	Middle Name	Office Name	Office Address	Email Address	Phone	Fax
Dr. Tuesday	Monday	H	Monday Office	1254 Humble [Berlin, DE]	woo1@lento-medical.com	4902495000	1

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Case List

Case #	Surgeon Name	Patient Name	DOPO	DOS	Surgery Type	Surgery Side	Status
K4fen6L	Tuesday, Monday	Sunday, Friday	Jul 27, 2020	Sep 25, 2020	TKR	Left	Ready to Upload MRI

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K4fen6L
Surgery Side: Left

MRI Upload Ready


Drop Zone
Click and Drag Patient MRI Files Here

Surgeon
First Name: Monday
Middle Name: H
Last Name: Tuesday
Office Name: Monday Office
Email Address: woo1@lento-medical.com
Phone: 4902495000

Patient
First Name: Friday
Middle Name: H
Last Name: Sunday
Sex: M
Age: 62
Date of Purchase Order: Jul 27, 2020
Date of Surgery: Sep 25, 2020
Surgery Type: TKR

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Surgeon: Office Name: Office Address: Email Address: Phone Number: Fax Number:	Seidel, Markus Dr. Markus Seidel Kastanienallee 2, Berlin, 10435, DE u@lentomedical.com 49304483568	NEW ORDER CURRENT ORDERS COMPLETED ORDERS VIEW PREFERENCES
---	---	---

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Current Orders

Case #	Patient Name	Order Date	Surgery Date	Surgery Type	Surgery Side	Status
K1gbQR	Park, Se Jung	Jul 27, 2020	Nov 17, 2020	TUR	Right	Surgeon Approval Ready

Case ID: K1gbQR
Status: Surgeon Approval Ready

COMPLETED

Patient

First Name: Se Jung
Last Name: Park
Date of Birth: Nov 27, 2009
Sex: Female
Age: 68
Date of Surgery: Nov 17, 2020
Insurance Information: STROVER HEALTHCO
Referral: CHARLSE UNIVERSITÄT MEDICIN
Referral Express: Refers express

Surgery Type: TUR
Surgery Side: Right
Default Implant: NEURONAL ALIGNMENT
Date of Purchase Order: Jul 27, 2020

Surgeon

Name: Seidel, Markus
Office Name: Dr. Markus Seidel
Address: Kastanienallee 2
Email Address: u@lentomedical.com
Office Number: 49304483568
Mobile Number: 49304483568
Fax: 49304483568

Distributor

Company Name: Refers express
Contact Person: Refers, Refers
Address: Kastanienallee 2
Email Address: u@lentomedical.com
Company Number: 49304483568
Mobile Number: 49304483568
Fax: 49304483568

Hospital

Hospital Name: Charité Universitätsmedizin
Contact Person: Refers, Refers
Address: Charité 1
Email Address: u@lentomedical.com
Hospital Number: 49304483568
Mobile Number: 49304483568
Fax: 49304483568

MRI Center

Center Name: Refers, Refers
Contact Person: Refers, Refers
Address: Kastanienallee 2
Email Address: u@lentomedical.com
Center Number: 49304483568
Mobile Number: 49304483568
Fax: 49304483568

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K1gbu1024L Total Knee Arthroplasty (TKA)

ZIMMER PERSONA

Ordering Surgeon: Dr. Surgeon, Lento
Patient Name: One, John
Date of Surgery: Nov 28, 2024

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

FEM IMP SIZE: 6
TIB IMP SIZE: E
Knee plate

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CAUTION!

This is to confirm that you are rejecting the pre-planned adjusted measurements made on the surgeon approval page. Please fill out rejection form below.

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